
# MODUL PRAKTIKUM

## KIMIA KLINIK 1

UNTUK KALANGAN SENDIRI

**LABORATORIUM PATOLOGI KLINIK  
FAKULTAS ILMU KESEHATAN  
UNIVERSITAS MUHAMMADIYAH SURABAYA 2019**

# MODUL PRAKTIKUM

## KIMIA KLINIK 1

**PENYUSUN :**

**KETUA : RAHMA WIDYASTUTI, S.Si, M.Kes**

**ANGGOTA : NUR VITA PURWANINGSIH, SST., M.Kes.**

**LABORATORIUM PATOLOGI KLINIK  
FAKULTAS ILMU KESEHATAN  
UNIVERSITAS MUHAMMADIYAH SURABAYA 2019**

## **VISI**

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**2.4 METHOD :** creatinine forms in alkaline solution an orange-red coloured complex with picric acid. The absorbance of this complex is proportional to the creatine concentration in the sample.

### 2.5 REACTION PRINCIPLE:

creatinine + picric acid  $\longrightarrow$  creatinine-picrate complex

---

## 2.6 REAGEN COMPOSITION :

|      |                             |                              |
|------|-----------------------------|------------------------------|
| PIC  | : 1x100 ml picric acid      | 26 mmol/l                    |
| NaOH | : 1x100 ml sodium hydroxide | 1,6 mol/l                    |
| STD  | : 1x25 ml standart          |                              |
|      | Creatinine                  | 2 mg/dl or 176,8 $\mu$ mol/l |

## 2.7 REAGENT PREPARATION

Dilute NaOH with dist. Water in the ratio 1+4 store the solution in a plastic bottle

Mix PIC and dilute NaOH for the working reagent in the ratio 1+1

The standart is ready for use

## 2.8 REAGENT STABILITY

The reagent / diluted sodium hydroxide are stable, even after opening, up to the stated expiry date when stored at 15...25  $^{\circ}$ C

Contamination must avoided

The working reagent protected from light is stable for 4 weeks at is 5...25  $^{\circ}$ C

## 2.9 SPECIMEN : serum, heparine plasma or urine

Avoid hemolysis!

Stability: 24 hours at 2-8  $^{\circ}$ C

Dilute urine : 1+49 with dist. Water

## 2.10 ASSAY:

Wavelength : Hg 492 nm (490 -510 nm)

Temperature : 25  $^{\circ}$ C (for 37  $^{\circ}$ C procedure ask Human GmbH)

Measurement : against air (increasing absorbance)

Warm the reagents and cuvettes, to the desired temperature and keep constant ( $\pm 0,5$  °C) for the duration of the test.

## 2.11 PROCEDURE

| Pipette cuvettes | Semi –micro  | Macro        |
|------------------|--------------|--------------|
| Sample/ standart | 100 $\mu$ l  | 200 $\mu$ l  |
| Working reagent  | 1000 $\mu$ l | 2000 $\mu$ l |

Mix and start the stopwatch, after 30 sec, read the absorbance A1 read the absorbance A2 exactly after 2 min  $A2-A1 =$  Absorbance sample or Absorbance standart

## 2.12 CALCULATION :

### 1. Serum/ plasma

Please use only the standart supplied with the kit.

$$C = 2,0 \times \frac{\Delta A \text{ sample}}{\Delta A \text{ standart}} \quad (\text{mg/dl})$$

$\Delta A$  standart

$$C = 176,8 \times \frac{\Delta A \text{ sample}}{\Delta A \text{ standart}} \quad (\mu\text{mol/l})$$

$\Delta A$  standart

### 2. Urine

$$C = 100 \times \frac{\Delta A \text{ sample}}{\Delta A \text{ standart}} \quad (\text{mg/dl})$$

$\Delta A$  standart

Creatinine concentration in 24 h urine:

$$C = \text{mg/dl} \times \text{ml urine} / 24 \text{ h} \times 0,01 \quad (\text{mg}/24\text{h})$$

$$C = \text{mg}/24 \text{ h} \times 0,00884 \quad (\text{mmol}/24 \text{ h})$$

$$\text{Creatinine} = \frac{\text{mg creatinine} / \text{dl urine ml urine} / 24 \text{ h}}{\text{mg creatinine} / \text{dl serum}} \times (\text{ml}/\text{min})$$

$$\text{mg creatinine} / \text{dl serum} \times 1440$$

conversion of (mg/dl) into ( $\mu\text{mol/l}$ ) and vice versa:

$$(\text{mg/dl}) \times 88,402 = (\mu\text{mol/l})$$

$$(\mu\text{mol/l}) \times 0,0113 = (\text{mg/dl})$$

### 2.13 Performance characteristic

linearity

the test is linear up to a creatinine concentration in serum of 13 mg/ dl or 1,150  $\mu\text{mol/l}$ , in urine of 500 mg/dl or 44,200  $\mu\text{mol/l}$ .

dilute samples with a higher concentration in serum, plasma or diluted urine 1+5 with physiological saline (0,9%) and repeat the assay. Multiple the result by 6.

Typical performance data can be found in the verification report accessible via.

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#### 3.4 METHOD :

urea is hydrolysed in the presence of water and urease to produce ammonia and carbon dioxide. In a modified berthelot reaction the ammonium ions react with hypochlorite and salicylate to form a green dye. The absorbance increase at 578 nm is proportional to the urea concentration in the sample.

#### 3.5 REAGEN :

|                                  |            |
|----------------------------------|------------|
| RGT 1: Phosphate buffer (pH 7,0) | 120 mmol/l |
| Sodium salicylate                | 60 mmol/l  |
| Sodium nitroprusside             | 5 mmol/l   |
| EDTA                             | 1 mmol/l   |

---

|  |                           |
|--|---------------------------|
| Rgt 2 : Phospate buffer (pH<13)  | 120 mmo/l                 |
| Hyphoclorite   | = 0,6 g/l Cl              |
| Irritates e yes and skin. Keep out of reach of children. Upon contact with the eyes, rinse thoroughly with water and consult a doctor. |                           |
| ENZ : enzyme   |                           |
| Urease   | >500 KU/l                 |
| STD : standart   |                           |
| Urea   | 80 mg/dl or 13,3mmol/l    |
| Equivalent to BUN  | 37,28 mg/dl or 6,2 mmol/l |
| Sodium azide   | 0,095%                    |

### 3.6 Reagent Prepration

RGT 2 and STD are ready for use

The enzyme reagent 1a is prepared by mixing the content of bottle ENZ with bottle

RGT1 : 1 ml ENZ + 100 ml RGT 1 or

10 l ENZ + 1000 ml RGT1

### 3.7 Reagent stability

The reagents are stable up to the stated expiry date when sealed and stored at 2...8 °C.

RGT 1, RGT 2, and ENZ are atable after opening for 6 weeks at 2...8 °C or 2 weeks at 15....25°C .

STD is stable up to the expiry date even after opening

The enzyme reagent 1a is a stable for 4 weeks at 2...8 °C or 2 weeks at 15....25°C .

Contamination after opening must be avoided

### 3.8 Specimen

Serum, plasma, except ammonium heparinate plasma and urine.

Dilute urine 1+100 with distilled water.

Do not use lipemic sera.

Do not use lipemic sera.

Serum or plasma can be stored for up to 3 days at 4<sup>0</sup>C, for longer periods the should be kept frozen at-20<sup>0</sup>C.

### 3.9 Assay

Wavelength : Hg 578 nm, 570- 600 nm

Optical path : 1 cm

Temperature : 20...25<sup>0</sup>C or 37<sup>0</sup>C

Measurement : against reagen blank. Only one reagent blank per series is required.

### 3.10 Procedure

| Pipette into cuvette   | Reagent<br>Blank | Sample or<br>Standart |
|--|------------------|-----------------------|
| Sample/ STD  | .....            | 10 µl                 |
| Enzyme reagent 1a  | 1000 µl          | 1000 µl               |
| Mix and incubate for 5 min at 20..25 <sup>0</sup> Cor for 3 min at 37 <sup>0</sup> C   |                  |                       |
| RGT 2  | 1000 µl          | 1000 µl               |
| Mix incubate for 10 min at20..25 <sup>0</sup> Cor for 5 min at 37 <sup>0</sup> C, measure the absorbance of the sample and the standart against the reagent blank within 60 min. |                  |                       |

### 3.11 Calculation of urea and BUN concentration

$$C = \frac{\Delta A \text{ sample}}{\Delta A \text{ standart}} \times \text{factor}$$

| Factor           | C(urea) |        | C(BUN) |        |
|------------------|---------|--------|--------|--------|
|                  | mg/dl   | mmol/l | mg/dl  | mmol/l |
| For serum/plasma | 80      | 13,3   | 37,28  | 6,2    |
| For urine        | g/l     | mmol/l | g/l    | mmol/l |
|                  | 80,8    | 1343   | 37,68  | 626,2  |

### 3.12 Conversion factor for BUN, urea

$$C \text{ (BUN)} = 0,466 \times C \text{ (urea)}$$

$$C \text{ (urea)} = 2,14 \times C \text{ (BUN)}$$

### 3.13 Performance characteristic

linearity

serum/plasma : up to 400 mg/dl or 66,6 mmol/l (urea)

urine : up to 400 g/l or 6600 mmol/l

samples with a higher urea concentration have to be diluted 1+1 with distilled water.

Repeat the assay and multiply the results by 2.

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

**4.4 METODE:** determination of uric acid by reaction with uricase. The formed H<sub>2</sub>O<sub>2</sub> reacts under catalysis of peroxidase with 3,5-dichloro-2-hydroxybenzenesulfonic (DCHBS) and 4-aminophenazone (PAP) to give a red-violet quinoneimine dye as indicator.

### 4.5 PRINSIP

### 4.6 KOMPOSISI REAGEN

|                           |            |
|---------------------------|------------|
| Buffer phospat (pH : 7,0) | 50 mmol/l  |
| 4-aminofenazon            | 0,3 mmol/l |
| DCHBS                     | 4 mmol/l   |
| Uricase                   | >200 U/l   |
| Peroxidase                | >1kU/l     |

#### 4.7 PROSEDUR

|   |         |               |
|---|---------|---------------|
| Pipette into cuvettes   | Rb      | Sample or std |
| Sample or std   | -----   | 20 µl         |
| Rgt   | 1000 µl | 1000 µl       |
| Mix incubate 10 min, at 20...25C or 5 min at 37C. Measure the absorbance of the sample / std against the reagent blank within 15 min. |         |               |

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

**5.4 METODE** : kinetic method for the determination of ASAT activity according to the recommendations of the expert panel of the IFCC (international federation of clinical chemistry) without pyridoxal phosphate activation.

### 5.5 Reaction principle

## 5.6 Reagent

|     |                          |            |
|-----|--------------------------|------------|
| BUF | : buffer/ enzyme reagent |            |
|     | TRIS buffer (pH 7,8)     | 100 mmol/l |
|     | L-aspartate              | 300 mmol/l |
|     | LDH                      | ≥ 0,9 kU/l |
|     | MDH                      | ≥ 0,6 kU/l |
| SUB | : substrate              |            |
|     | 2-oxoglutarate           | 60 mmol/l  |
|     | NADH                     | 0,9 mmol/l |

## 5.7 Reagent preparation and stability

REF : 12211 : pipette 1 ml from bottle SUB into one bottle BUF mix thrughly  
REF : 12011 : pipette 2 ml from bottle SUB into one bottle BUF mix thrughly  
The working reagent is stable for 4 weeks at 2...8 °C and 5 days at 15....25 °C

## 5.8 Specimen

Serum, heparinised plasma or EDTA plasma  
Avoid hemolysys !

## 5.9 Assay

Wavelength : Hg 365 nm, 340 nm or Hg 334 nm  
Optical path : 1 cm  
Temperature : 25 °C, 30°C, or 37 °C  
Measurement : against air (decreasing absorbance)  
Warm the reagents and the cuvettes to the desired temperature must be constant (±0,5°C) for the duration of the test.

## 5.10 Procedure :

| Pipette info cuvettes  | 25°C, 30°C | 37°C    |
|--|------------|---------|
| Sample   | 200 µl     | 100 µl  |
| Working reagent  | 1000 µl    | 1000 µl |
| Mix, read the absorbance after 1 minute and at the same time start the stop watch.<br>Read the absorbance again exactly after 1,2 and 3 minutes. |            |         |

### 5.11 Calculation

For A/min within 0,06-0,08 (Hg 365 nm) or 0,12-0,16 (Hg 334 nm, 340 nm) (procedure 1+2) use only measurements from the first 2 minutes for calculation (1 min incubation, 2 min measurement)

| U/l = A/min x | Sample start |      | Reagent start |      |
|---------------|--------------|------|---------------|------|
|               | 25°C, 30°C   | 37°C | 25°C, 30°C    | 37°C |
| Hg 334 nm     | 971          | 1780 | 1173          | 2184 |
| 340 nm        | 952          | 1745 | 1151          | 2143 |
| Hg 365 nm     | 1765         | 3235 | 2132          | 3971 |

Conversion factor from traditional units (U/l) in SI-units (kat/l):

1 U/l =  $16,67 \times 10^3$   $\mu$ kat/l

1  $\mu$ kat/l = 60 U/l

### 5.12 Performance characteristic

Linearity

If the absorbance change per minute (A/ min) or the activity exceed

| Wavelength (nm) | A/min | 25°C, 30°C (U/l) | 37°C (U/l) |
|-----------------|-------|------------------|------------|
| Hg 365          | 0,080 | 170              | 320        |
| Hg 334/340      | 0,160 | 190              | 350        |

Dilute 0,1 ml of the sample with 0,9 ml physiological saline (0,9%) and repeat the assay using this dilution. Multiply the result by 10.

In sera with very high activities, the initial absorbance may be very low as most of the NADH may have been consumed before the first reading. In this case rerun the sample after dilution as described above.

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**6.4 METODE** : kinetic method for the determination of ALAT activity according to the recommendations of the expert panel of the IFCC (international federation of clinical chemistry) without pyridoxalphosphate activation.

### 6.5 REACTION PRINCIPLE:

2 Oxoglutarate + L-alanine  $\xrightarrow{\text{GPT}}$

GPT

L-glutamate + pyruvate

Pyruvate + NADH + H<sup>+</sup> +

### 6.6 KOMPOSISI REAGEN

BUF = buffer / enzyme reagent

TRIS buffer (pH: 7,5)

150 mmol/l

L-alanine

750 mmol/l

LDH

≥1,2 kU/l

---

SUB = Substrate  
2-oxoglutarate 90 mmol/l  
NADH 0,9 mmol/l

### 6.7 PERSIAPAN REAGEN

REF 1202 : pipette 2 ml from bottle SUB into one bottle BUF mix thoroughly

### 6.8 PROSEDUR

| Pipette into cuvettes | 25C     | 37C     |
|-----------------------|---------|---------|
| Sample                | 200 µl  | 100 µl  |
| Working reagent       | 1000 µl | 1000 µl |

Mix, read the absorbance after 1 minute and at the same time start the stop watch.  
Read the absorbance again exactly 1,2 and 3 minute

### 6.9 HARGA NORMAL

| Temperature | 25 C   | 30C    | 37C |
|-------------|--------|--------|-----|
| Men up to   | 22 U/l | 30 U/l | 42C |
| Women up to | 17 U/l | 23 U/l | 32C |

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
### 7.4 PRINSIP

Bilirubin reacts with sulphanic acid (DSA) to form a red azo dye. The absorbance of this dye at 546 nm is directly proportional to the bilirubin concentration in the presence of an accelerator : total bilirubin = direct + indirect bilirubin.

### 7.5 CALCULATION

Calculate the concentration of total and direct bilirubin by using the factor 13,0

Bilirubin concentration (mg/dl) = absorbance x 13,0

(mg/dl) x 17,1 = (µmol/l)

## 7.6 Normal values

| Total bilirubin         | (mg/dl) |
|-------------------------|---------|
| At birth up to          | 5       |
| 5 days up to            | 12      |
| 1 month up to           | 1,5     |
| Adults up to            | 1,1     |
| <b>Direct bilirubin</b> |         |
| Adults up to            | 0,25    |

## 7.7 PEMERIKSAAN BILIRUBIN TOTAL PROSEDUR

|  |              |              |
|--|--------------|--------------|
| Pipette  | Sample blank | Sample       |
| TBR  | 1000 $\mu$ l | 1000 $\mu$ l |
| TNR  | -----        | 1 drop       |
| Mix thoroughly incubate for 5 min  |              |              |
| Sample   | 100 $\mu$ l  | 100 $\mu$ l  |
| Mix incubate at room temperature for 10 to 30 min. measure the absorbance of sample against sample blank |              |              |

1 drop  $\approx$  40  $\mu$ l

## 7.8 PEMERIKSAAN BILIRUBIN DIRECT PROSEDUR

|   |              |              |
|---|--------------|--------------|
| Pipette   | Sample blank | Sample       |
| DBR   | 1000 $\mu$ l | 1000 $\mu$ l |
| DNR   | -----        | 1 drop       |
| Mix thoroughly incubate for 2 min   |              |              |
| Sample  | 100 $\mu$ l  | 100 $\mu$ l  |
| Mix incubate at room temperature for exactly 5 min. measure the absorbance of sample against sample blank |              |              |

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**8.4 METODE :** cupric ions react with protein in alkaline solution to form a purple complex. The absorbance of this complex is a proportional to the protein concentration in the sample.

### 8.5 KOMPOSISI REAGEN

|                           |            |
|---------------------------|------------|
| RGT : sodium hydroxide    | 200 mmol/l |
| Potassium sodium tartrate | 32 mmol/l  |
| Copper sulfate            | 12 mmol/l  |
| Potassium iodide          | 30 mmol/l  |
| Irritant R 36/38          |            |
| STD                       |            |
| Protein                   | 8 g/dl     |
| Or                        |            |
| Sodium azide              | 0,095%     |

### 8.9 PROSEDUR

| Pipette into cuvettes | Reagent Blank | Sampe/standart |
|-----------------------|---------------|----------------|
| Sampe/standart        | -----         | 20 µl          |
| Reagent               | 1000 µl       | 1000 µl        |

Mix incubate for 10 min at 20-25 C. measure the absorbance of the sample and standart against the reagent blank within 30 min.

---

### 8.10 Calculation

1. With factor :  
 $C = 19 \times \text{absorbance (g/dl)}$
2. With standart  
 $C = 8 \times \frac{\text{absorbance sample}}{\text{Absorbance standart}} \text{ (g/dl)}$

### 8.11 TAMPILAN KARAKTERISTIK

The tes is linear up to a protein concentration of 12 g/dl or 120 g/l. dilute samples with a higher concentration 1+1 with physiological saline (0,9%). Multiply the result by 2.

### 8.12 Normal range in serum and plasma

|                                  |              |
|----------------------------------|--------------|
| Normal born babies               | 4,6-7,0 g/dl |
| Children from 3 years and adults | 6,6-8,7 g/dl |

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### 9.3 METODE

Bromocresol green with albumin in citrate buffer a coloured complex. The absorbance of this complex is proportional to the albumin concentration in the sample.

### 9.4 PROSEDUR

|  |               |               |
|--|---------------|---------------|
| Pipette into cuvettes  | Reagent blank | Sample or std |
| Sample or std  | -----         | 10 µl         |
| Rgt  | 1000 µl       | 1000 µl       |
| Mix incubate for 5 min, at 20..25C. measure the absorbance of the sample and the standart against the reagent blank within 30 min. |               |               |

### 9.5 PERHITUNGAN

$$C = 4 \times \frac{\text{Absorbance sample (g/dl)}}{\text{Absorbance std}}$$

### 9.6 TAMPILAN KARAKTERISTIK

The tes is linear up to an albumin concentration of 7 gr/dl. Dilute sample with a higher concentration 1+1 with physiological saline (0,9%) multiply the result by 2.

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### 10.4 METODE

“Optimized standart method” according to the recommendation of the german clinical chemnistry association.

## 10.5 REACTION PRINSIP

## 10.6 KOMPOSISI REAGEN

|     |                                       |  |              |
|-----|---------------------------------------|--|--------------|
| BUF | Buffer                                |  |              |
|     | Diethanolamine buffer (pH 10,35 ±0.2) |  | 1.25 mol/l   |
|     | Magnesium chloride                    |  | 0.625 mmol/l |
| SUB | Substrate                             |  |              |
|     | p-Nitrophenyl phosphate               |  | 55 mmol/l    |

## 10.7 PERSIAPAN REAGEN

REF 12017 : pipette 2 ml from bottle SUB into one bottle BUF mix throughly  
The working reagent is stable for 4 weeks at 2...8C, for 5 days at 15...25  
C. the working reagent must be kept light protected.

## 10.8 PROSEDUR

|   |                     |
|---|---------------------|
| <b>Pipette into cuvette</b>   | <b>25C, 30C,37C</b> |
| Sample  | 20µl                |
| Working reagent   | 1000 µl             |
| Mix read the absorbance after 1 min and at the same time start the absorbance. Read the absorbance again exactly after 1,2 and 3 min. |                     |

## 10.9 TAMPILAN KARAKTERISTIK

If the absorbance change per minute (Abs/min) exceeds 0.250 dilute 0.1 ml of the sample with 0.5 ml physiological saline (0,9%) and repeat the assay using this dilution. Multiply the result by 6.

---

## 10.10 HARGA NORMAL

| Temperature                | 25C<br>U/l | 30C<br>U/l | 37C<br>U/l |
|----------------------------|------------|------------|------------|
| Women                      | 40-190     | 49-232     | 63-306     |
| Men                        | 50-190     | 61-232     | 80-306     |
| Children up to 15<br>years | Up to 400  | Up to 488  | Up to 644  |
| Children up to 17<br>years | Up to 300  | Up to 366  | Up to 483  |

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### 11.4 METODE

The trigliserides are determined after enzymatic hydrolysis with lipases. Indicator is quinoneimine formed from hydrogen peroxide, 4amino antipyrine and 4 chlorophenol under the catalytic influence of peroxidase

### 11.5 REACTION PRINCIPLE

Triglycerides  $\longrightarrow$  glycerol + fatty acids

Glycerol + ATP  $\longrightarrow$  glycerol-3-phosphate + ADP

glycerol-3-phosphate + O<sub>2</sub>  $\longrightarrow$  dihydroxyacetone phosphate + H<sub>2</sub>O<sub>2</sub>

H<sub>2</sub>O<sub>2</sub> + 4 aminoantipyrine  $\longrightarrow$  quinoneimine + HCL + H<sub>2</sub>O + 4-chlorophenol

### 11.6 KOMPOSISI REAGEN

|                             |             |
|-----------------------------|-------------|
| RGT : Pipes buffer (pH 7,5) | 50 mmol/l   |
| 4-cholophenol               | 5 mmol/l    |
| 4-aminoantipyrine           | 0,25 mmol/l |
| Magnesium ions              | 4,5 mmol/l  |
| ATP                         | 2 mmol/l    |
| Lipases                     | ≥1,3 U/ml   |
| Peroxidase                  | ≥0,5 U/ml   |
| Glycerol kinase             | ≥ 0,4 U/ml  |
| Glycerol-3-phospate oxidase | ≥ 1,5 U/ml  |
| STD                         |             |
| Triglyserides               | 200 mg/dl   |

### 11.7 PROSEDUR

| Pipette into cuvettes   | RB      | Sample or STD |
|---|---------|---------------|
| Sample/STD  | -       | 10 µl         |
| RGT   | 1000 µl | 1000 µl       |
| Mix and incubate for 10 minute at 20...25C or for 5 minute at 37C. measure the absorbance of the sample and the standart against the reagent blank within 60 min. |         |               |

### 11.8 PERHITUNGAN

$$200 \times \frac{\text{Abs sample}}{\text{Abs standart}} \text{ (mg/dl)}$$

### 11.20 TAMPILAN KARAKTERISTIK

The test is linear up to a trigliserides concentration of 1000 mg/dl Or 11,4 mmol/l. sample with a higher contraction have to be diluted 1+4 with phsiogical saline (0,9%) and retested. Multiply the result by 5.

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### 12.4 METODE

The cholesterol is determined after enzymatic hydrolysis and oxidation. The indicator quinoneimine is formed from hydrogen peroxide and 4-aminophenazone in the presence of phenol and peroxidase.

### 12.5 REACTION PRINCIPLE

Cholesterolester + H<sub>2</sub>O<sub>2</sub> → kolesterol + fatty acid

Cholesterol + O<sub>2</sub> → cholestene -3-one + H<sub>2</sub>O<sub>2</sub>

2H<sub>2</sub>O<sub>2</sub> + 4aminophenazone + phenol → quinoneimine + 4 H<sub>2</sub>O

## 12.6 KOMPOSISI REAGEN

|                                  |            |
|----------------------------------|------------|
| RGT : Phosphate buffer (pH: 6,5) | 100 mmol/l |
| 4-aminophenazone                 | 0,3 mmol/l |
| Phenol                           | 5 mmol/l   |
| Peroxidase                       | >5 KU/l    |
| Cholesterolesterase              | >150U/l    |
| Cholesteroloxidase               | >100 U/l   |
| Sodium azide                     | 0,05%      |
| STD                              |            |
| Cholesterol                      | 200 mg/dl  |

## 12.7 PROSEDUR

| Pipette into cuvettes  | Reagent blank | Sample or standart |
|--|---------------|--------------------|
| Sample/standart  | -             | 10 $\mu$ l         |
| reagent  | 1000 $\mu$ l  | 1000 $\mu$ l       |
| Mix incubate 10 min at 20...25C or 5 min at 37C. measure the absorbance of the sample / std against the reagent blank within 60 min. |               |                    |

## 12.8 PERHITUNGAN

$$C = \frac{200 \times \text{Abs sample}}{\text{Abs standart}} \quad (\text{mg/dl})$$

## 12.9 TAMPILAN KARAKTERISTIK

The test is linear up to a cholesterol concentration of 750 mg/dl. Dilute samples with a higher cholesterol concentration 1+2 with physiological saline (0,9%) and repeat

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

### 13.4 METODE

The glucose is determined after enzymatic oxidation in the presence of glucose oxidase. The formed hydrogen peroxidase reacts under catalysis of peroxidase with phenol and 4-aminophenazone to a red violet quinoneimine dye as indicator.

### 13.5 REACTION PRINCIPLE

### 13.6 KOMPOSISI REAGEN

|            |                            |             |
|------------|----------------------------|-------------|
| <b>RGT</b> | phosphate buffer (pH: 7,5) | 0.1 mol/l   |
|            | 4-aminophenazone           | 0.25 mmol/l |
|            | Phenol                     | 0.75 mmol/l |
|            | glucose oxidase            | >15KU/l     |
|            | Peroxidase                 | >1.5 KU/l   |
|            | Mutarotase                 | >2.0 KU/l   |
|            | Stabilizer                 |             |
| <b>STD</b> | glucose                    | 100 mg/dl   |

### 13.7 PROSEDUR

|                      | Semi micro    |               |
|----------------------|---------------|---------------|
| Pipette into cuvette | STD or sample | Reagent blank |
| STD or sample        | 10 $\mu$ l    | -             |

---

|  |        |        |
|--|--------|--------|
| RGT  | 1000µl | 1000µl |
| Mix incubate for 10 min at 20...25C or 5 min at 37C. measure the absorbance of STD and the sample against the reagent blank within 60 min. |        |        |

### 13.8 PERHITUNGAN

$$C = 100 \times \frac{\text{Abs sample (mg/dl)}}{\text{Abs standart}}$$

### 13.9 TAMPILAN KARAKTERISTIK

The test is linear up to a glucose concentration of 400 mg/dl. Dilute the sample 1+2 with dest water, if the glucose concentration of the sample is over this limit and repeat the determination. Multiply the rest by 3.

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