MODUL PRAKTIKUM

KIMIA KLINIK 2

UNTUK KALANGAN SENDIRI

LABORATORIUM PATOLOGI KLINIK
FAKULTAS ILMU KESEHATAN
UNIVERSITAS MUHAMMADIYAH SURABAYA 2018

MODUL PRAKTIKUM

KIMIA KLINIK 2

PENYUSUN

KETUA : RAHMA WIDYASTUTI, S.Si, M.Kes

ANGGOTA : NUR VITA PURWANINGSIH, S.ST,M.Kes

**LABORATORIUM PATOLOGI KLINIK
FAKULTAS ILMU KESEHATAN
UNIVERSITAS MUHAMMADIYAH SURABAYA 2018**

VISI

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Fitrotin Aizah, S.ST, M.Si

Surabaya, September 2018
PJMK

A handwritten signature in blue ink, which appears to read 'Rahma Widayastuti'. It includes a double quote mark above the signature and a horizontal line underneath.

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METHOD : Precipitant

PRINCIP : The chylomicrons, VLDL (Very low density lipoproteins) and LDL (low density lipoproteins) are precipitated by addition of phosphotungstic acid and magnesium chloride. After centrifugation the supernatant fluid contains the HDL (High Density Lipoprotein) fraction , which is assayed for HDL Cholesterol with the human cholesterol liquicolor test.

SPECIMEN : Serum, Plasma Heparin atau EDTA

REAGEN COMPOSITION

PREC : phosphotungstic acid
magnesium chloride

STD : Standard Cholesterol

REAGENT PREPARATION :

Macro : use undiluted (Reagen PREC Ca)

SemiMicro : Diluted the contents of 1 bottle PREC with 20 ml Distilled water, or dilute 4 parts of the bottle with 1 part distilled water (4 +1)

PROCEDURE:

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1. Precipitasi

Pipette into centrifuge tubes	Macro	Semi-Micro
Sample	500 µl	200 µl
Prec Ca	1000 µl	---
Prec Cb	---	500µl

Mix well, incubate for 10 minute at room temperature. Centrifuge for at least 2 minutes at 10000 g, alternatively for 10 minutes at 4000 g.

After centrifugation separate the clear supernatant from the precipitate within 1 hours and determine the cholesterol concentration using cholesterol reagent.

2. Cholesterol Determination

Pipette into cuvettes	Reagent blank	STD	Sample
Dist water	100 µl	---	---
STD	---	100 µl	---
HDL Supernatant	---	---	100 µl
Reagent	1000 µl	1000 µl	1000 µl

Mix incubate for 5 minutes at 37 C or 10 minutes at 20...25C. Measure the absorbance of the sample and the standart, respeively, against the reagent balnk within 60 minutes

CALCULATION :

1. Macro Method

$$C = \frac{150 \times \text{Abs Sample mg/dl}}{\text{Abs STD}}$$

2. Semi micro method

$$C = \frac{175 \times \text{Abs Sample mg/dl}}{\text{Abs STD}}$$

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The LDL cholesterol concentration (LDL-C) is calculated from the total cholesterol concentration (TC-C), the HDL Cholesterol contration (HDL-C) and the triglycerides concentration (TG) according to Friedewald.

Calculation

$$\text{LDL-C} = \text{TC} - \text{HDL-C} - \frac{\text{TG}}{5}$$

Clinical Interpretation

1. Suspicious : 150 mg/dl
2. Elevated : 190 mg/dl

Reference :

Gordon, T and M. amer. J.Med. 62 707 (1977)

Friedewald, W.T.et.al. Clin.Chem. 18,499 (1972)

ISO 15223 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied

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METHOD

Sodium is precipitated with Mg-uranl acetate ; the uranyl ions remaining in suspension form a yellow – brown complex with thio glcolic acid. The difference between reagent blank (without precipitation of sodium) and anlsis is proportional to the sodium concentration.

REAGEN COMPOSITION :

For 20 macro / 60 semimicro determinations

<u>Prec</u>	60 ml prepipitating solution	
	Uranyl acetate	19 mmol/ l
	Magnesium acetate	140 mmol/l
RGT	60 ml colour reagent	
	Ammonium thioglycolate	550 mmol/ l
	Ammonia	550 mmol/ l
STD	2 ml standart	
	Sodium (Na+)	150 mmol/l

SPECIMEN : serum

PROCEDURE

(ul)	Makro			Semi – micro		
	RB	STD	SPL	RB	STD	SPL
R4/Std	-	20	-	-	10	-
Sample	-	-	20	-	-	10
R1	-	2000	2000	-	1000	1000
Close tubes and mix well. Allow to stand for 5 minutes. Shake intensively for ca. 30 second. Allow to stand for 30 minutes. Centrifuge at high speed for 5 – 10 minutes.						
Supernatant	-	50	50	-	20	20
R1	50	-	-	20	-	-
R2	2000	2000	2000	1000	1000	1000
Mix well. After 10 -15 minutes, measure Abs. Blank (dARB), Standart (Dastd) And sample (Daspl) Against Water At 360 – 410nm (Hg 365 or Hg 405)						

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METHOD :

Pottassium ions in a protein – free alkaline medium react with sodium tetraphenylboron to produce a finely dispersed turbid suspension of pottassium tetraphenylboron. The turbidity produced is proportional to the pottassium concentration and read photometrically.

REAGENT PREPARATION

Mix the contant of bottle TPB with the contant of bottle NAOH. For smaller amounts of working reagent mix TPB and NAOH in a ratio of 1 + 1 (see reagent stability)

Allow to stand for 15-30 min prior to use

PREC and std are ready for use

STD is used undiluted directly in the determination

SPECIMEN : serum and lithium heparin plasma

ASSAY:

Wavelength : 578 nm, Hg 578 nm

Optical path : 1 cm

Temperature : 20 -25 C

Measurement : against reagent blank

Only one reagent blank per series is required

PROCEDURE

Precipitation				
Pipette into centrifuge tubes				
	Macro	Semi micro		
Specimen PREC	100 µl 1000 µl		50 µl 500 µl	
mix carefully, centrifuge at high speed for 5- 10 min				
Determination				
Pipette into cuvettes				
	STD	sample	STD	Sample
Working reagent STD Supernatant	2000 µl 200 µl ----	2000 µl ----- 200 µl	1000 µl 100 µl ---	1000 µl --- 100 µl
To produce a homogeneous turbidity the STD or the clear supernatant have to be added to the centre of the surface of the working reagent in the cuvette. Mix each cuvette carefully before proceeding to the next sample. Allow to stand at least for 5 min.				
Measure absorbance of the standart (STD) and the sample (SPLE) against working reagent blank between 5 and 30 minutes				

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Diagnostic reagent for quantitative in vitro determination of calcium in serum, plasma or urine on photometric systems

Order information

Cat. No.	Kit	Size	
1 1130 99 83 021	R	5x	25 ml +1x3 ml standart
1 1130 99 83 026	R	6x	100 ml
1 1130 99 83 023	R	1x	1000 ml
1 1130 99 83 704	R	8X	50 ml
1 1130 99 83 717	R	6X	100 ml
1 1130 99 83 030		6X	3 ml standart

Calcium plays an essential role in many cell functions : intracellularly in muscle contraction and glycogen metabolism, extracellularly, in bone mineralization, in blood coagulation and in transmission of nerve impulses. Calcium is present in plasma in three forms : free, bound to proteins or complexed with anions as phosphate, citrate and bicarbonate. Decreased total calcium levels can be associated with diseases of the bone apparatus (especially osteoporosis), kidney diseases (especially under dialysis), defective intestinal absorption and hypoparathyroidism. Increased total calcium can be measured in hypoparathyroidism, malignant diseases with metastases and sarcoidosis. Calcium measurements also help in monitoring of calcium supplementation mainly in the prevention of osteoporosis.

Method

Photometric test using arsenazo

Principle

Calcium with arsenazo III at neutral pH yields a blue colored complex, whose intensity is proportional to the calcium concentration. Interference by magnesium is eliminated by addition of 8-hydroxyquinoline-5-sulfonic acid

Reagents

Components and concentration in the test

Reagent :

Phosphate buffer pH 7,5 50 mmol/l

8-hydroxyquinoline-5-sulfonic acid 5 mmol/l

Arsenazo III 120 umol/l

Detergents

Standart : 10 mg/dl (2.5 mmol/l)

Storage instructions and reagents stability

The reagent and the standart are stable up to the end of the indicated month of expiry, stored at 2-8 °C

Warning and precautions

1. As calcium is an ubiquitous ion essential precaution must be taken against accidental contamination. Only use disposable materials.
2. Traces of chelating agent, such as EDTA can prevent the formation of the colored complex.
3. The reagent contains sodium azide (0,95 g/l) as preservative. Do not swallow ! avoid contact with skin and mucous membranes.
4. Take the necessary precautions for the use of laboratory reagents.

#Waste management

Please refer to local legal requirements

#Reagent preparation

The reagent and the standart are ready-to-use

#Materials required but not provided

NaCl solution 9 g/l

General laboratory equipment

#specimen

Serum, heparin plasma or urine

Do not use EDTA plasma.

Stability in serum/plasma :

7 days	at	20-25 °C
3 weeks	at	4-8 °C
8 months	at	-20 °C

Stability in urine :

2 days	at	20-25 °C
4 days	at	4-8 °C
3 weeks	at	-20 °C

Add 10 ml of concentrated HCl to 24 h urine and heat the specimen to dissolve calcium oxalate.

Discard contaminated specimens.

#Assay procedure

Application sheets for automated systems are available on request.

Wavelength 650 nm, Hg 623 nm (639-670 nm)

Optical path 1 cm

Temperature 20-25 °C / 37 °C

Measurement againts reagent blank

	Blank	Sample or standard
Sample or standard	-	10 ul
Dist. Water	10 ul	-
Reagent	1000 ul	1000 ul

Mix, incubate for 5 min, and read absorbance againts reagent blank

Calculation

With standart or calibrator

A sample

Calcum [mg/dl] = ----- x conc. Std / cal [mg/dl]

A std / cal

#conversion factor

Calcium [mg/dl] x 0.2495 = calcium [mmol/l]

#calibrators and controls

For the calibration of automated photometric systems the TruCal U calibrator is recommended. For internal quality control TruLab N and P control should be assayed with each batch of samples.

#Performance Characteristic

Measuring range

The test has been developed to determine calcium concentrations within a measuring range from 0.04 -25 mg/dl (0.01-6.25 mmol/l). When values exceed this range samples should be diluted 1+1 with NaCl solution (9 g/l) and the result multiplied by 2.

Specificity / interferences

No interference was observed by ascorbic acid up to 30 mg/dl, bilirubin up to 40 mg/dl, hemoglobin up to 500 mg/dl, lipemia up to 2,000 mg/dl triglycerides and magnesium up to 15 mg/dl.

Sensitivity

The lower limit of detection is 0.04 mg/dl (0.01 mmol/l)

#precision (at 20-25 °C)

Intra-assay precision N = 20	Mean [mg/dl]	SD [mg/dl]	CV[%]
Sample 1	8.79	0.09	1.04
Sample 2	12.5	0.15	1.20
Sample 3	14.0	0.24	1.73

Intra-assay precision N = 20	Mean [mg/dl]	SD [mg/dl]	CV[%]
Sample 1	8.82	0.18	2.01
Sample 2	12.3	0.11	0.90
Sample 3	13.7	0.26	1.92

method comparison

A comparison between calcium FS (y) and a commercially available test (x) using 70 samples gave following result : $y = 1.02 x - 0.20$; $r = 0.999$.

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RGT and STD are ready to use

REAGENT STABILITY

RGT and STD are stable even after opening up to the stated expiry date when stored at 2....25°C in the dark. Contamination must be avoided.

SPECIMEN

Serum, CSF, urine

Stability in serum at 2.....25°C 7 days

Stability in urine at 2.....8 °C 7 days

NOTES

Blood should be immediately processed otherwise an electrolyte exchange between serum and erythrocytes may lead to a false low chloride determination.

Precipitates in urine may be redissolved by warming up

ASSAY

Wavelength : 590 nm (590 -600 nm), Hg 578 nm

- Optical path : 1 cm
Temperature : 20....37 °C
Measurement : against reagent blank Only one blank per series is required

PROCEDURE

Macro method: predilute STD and the sample 1 + 40 with distilled water, e.g. 50 µl STD / sample + 2000 µl dist water

Semi- micro method: predilut STD and the sample 1 + 50 with distilled water, eg 20 µl STD / sample + 1000 µl dist water.

Use the preiluted STD/ sample for the following pipetting scheme

Pipette into cuvettes	Macro		Semi-micro	
	STD	Sample	STD	sample
STD	50 µl	----	20 µl	----
Sample	----	50 µl	----	20 µl
RGT	2000 µl	2000 µl	1000 µl	1000 µl

Mix, incubate for 5 minutes in the dark and measure the absorbance of sample (A sample) and STD (Astd) within 60 minutes against the reagent blank. Do not expose to light!

CALCULATION

$$C = 100 \times \underline{A \text{ sample}} \text{ (mmol/l)}$$

Astandart

$$C= 355 \times \underline{A \text{ sample}} \text{ (mg/dl)}$$

Astandart

The result may be converted according to the following equations:

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
METHOD : creatinine forms in alkaline solution an orange-red coloured complex with picric acid. The absorbance of this complex is proportional to the creatine concentration in the sample.

REACTION PRINCIPLE:

REAGEN COMPOSITION :

PIC : 1x100 ml picric acid 26 mmol/l

NaOH : 1x100 ml sodium hydroxide 1,6 mol/l

STD : 1x25 ml standart

Creatinine 2 mg/dl or 176,8 µmol/l

REAGENT PREPARATION

Dilute NaOH with dist. Water in the ratio 1+4 store the solution in a plastic bottle

Mix PIC and dilute NaOH for the working reagent in the ratio 1+1

The standart is ready for use

REAGENT STABILITY

The reagent / diluted sodium hydroxide are stable, even after opening, up to the stated expiry date when stored at 15...25 °C

Contamination must avoided

The working reagent protected from light is stable for 4 weeks at is 5...25 °C

SPECIMEN : serum, heparine plasma or urine

Avoid hemolysis!

Stability: 24 hours at 2-8 °C

Dilute urine : 1+49 with dist. Water

ASSAY:

Wavelength : Hg 492 nm (490 -510 nm)

Temperature : 25 °C (for 37 °C procedure ask Human GmbH

Measurement : against air (increasing absorbance)

Warm the reagents and cuvettes, to the desired temperature and keepnconstant ($\pm 0,5$ °C) for the duration of the test.

PROCEDURE

Pipette cuvettes	Semi –micro	Macro
Sample/ standart	100 µl	200 µl
Working reagent	1000µl	2000 µl
Mix and start the st opwach, after 30 sec, read the absorbance A1 read the absorbance A2 exatly after 2 min A2-A1 = Absorbance sample or Absorbance standart		

CALCULATION :

1. Serum/ plasma

Please use only the standart supplied with the kit.

$$C = 2,0 \times \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standart}}} \text{ (mg/dl)}$$

$\Delta A_{\text{standart}}$

$$C = 176,8 \times \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standart}}} \text{ (\mu mol/l)}$$

$\Delta A_{\text{standart}}$

2. Urine

$$C = 100 \times \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standart}}} \text{ (mg/dl)}$$

$\Delta A_{\text{standart}}$

Creatinine concertration in 24 h urine:

$$C = \text{mg/dl} \times \text{ml urine} / 24 \text{ h} \times 0,01 \text{ (mg/24h)}$$

$$C = \text{mg}/24 \text{ h} \times 0,00884 \text{ (mmol/24 h)}$$

$$\text{Creatinine} = \frac{\text{mg creatinine}}{\text{dl urine}} \times \frac{\text{ml urine}}{24 \text{ h}} \times \frac{(\text{ml/min})}{\text{mg creatinine} / \text{dl serum} \times 1440}$$

conversion of (mg/dl) into ($\mu\text{mol/l}$) and vice versa:

$$(\text{mg/dl}) \times 88,402 = (\mu\text{mol/l})$$

$$(\mu\text{mol/l}) \times 0,0113 = (\text{mg/dl})$$

Performance characteristic

linearity

the test is linear up to a creatinine concertration in serum of 13 mg/ dl or 1,150 $\mu\text{mol/l}$, in urine of 500 mg/dl or 44,200 $\mu\text{mol/l}$.

dilute samples with a higher concentration in serum, plasma or diluted urine 1+5 with physiological saline (0,9%) and repeat the assay. Multiple the result by 6.


Typical performance date can be found in the verification report accessible via.

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